## PRODUCT SPECIFICATION





### 1 GENERAL PRODUCT DESCRIPTION

#### Description

A personal eyewash bottle containing sterile, buffered, saline eye wash solution.

The Cederroth Eye Wash is intended for washing away foreign particles and dangerous liquids, like alkali or acids, from an eye in an emergency first aid situation at work places.

The Cederroth buffered Eye Wash also has a neutralizing effect on splashes from alkalis and acids. The effect is in general stronger for alkalis than for acids.

The Cederroth Eye Wash is intended for the first acute phase at the very accident site before getting to a plumbed or self-contained eye wash equipment or seeking medical attention.

The bottle is easily opened and ready to use by twisting off the dustcover. The integrated eye cup helps keeping the eye open during rinsing and ensures a generous flow without unnecessary waste.

Accessory: Wall bracket for Cederroth Eye Wash, REF 7200

# **Storage information**

Storage in +2°C - +40°C

### 2 DETAILED PRODUCT INFORMATION

## Content per consumer unit

<u>Product/Material</u>	<u>Qty</u>	Dimensions (b	ottle)
Cederroth Eye Wash	2	Diameter	6,6 cm
		Height	23,5 cm
		Fill volume	500 ml

#### **Material**

Liquid 0.9% saline solution, borate buffered to neutral pH. Preservative free.

Bottle Leak proof polypropylene bottle produced with Blow-Fill-Seal technique.

Eye cup in transparent polypropylene. Dust cover in green hard plastics (ABS).

**Shelf life:** 4,5 years

## **3 LOGISTIC INFORMATION**

	CU	SKU	Pallet
GTIN:	7310617252000	7310807262001	7312027252007
Content:	2 bottles	6 CU	11 SKU/layer x 5 layers
Size (w x l x h mm):	147 x 75 x 246	235 x 300 x 260	Pallet 800 x 1200 x 1450
Gross weight:	1196 g	8000 g	465 kg
Packing material:	Corrugated board 84g	Corrugated board 824 g	
	Plastic 112 g		



Version: 03 Issued: 2017-01-04